

# siRNA-27

## Application Guide

### Table of Contents

Package Contents and Storage Conditions .....	2
Storage Conditions .....	2
Related products .....	2
Additional materials recommended .....	2
Notice to purchaser .....	2
RNAi Collection Overview .....	3
Product Application .....	3
Reconstitution of siRNA duplexes in duplex buffer .....	3
Outline of an RNAi Experiment .....	4
Choice of Transfection Method .....	5
Transfection Optimization and protocol .....	6
Protocol for transient transfection in a 96 well plate .....	6
Positive and Negative Controls .....	7
FAQ: siRNA-27 kit .....	7
References .....	9



#### UK & Rest of World

184 Millon Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

#### Switzerland

Centra Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

#### Deutschland

Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 779099  
Fax: +49 (0) 69 13376880

#### North America

23591 El Toro Rd, Suite #180  
Lake Forest, CA 92630  
Tel: +1 800 987 0985  
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com  
AMS Biotechnology

## ***Package Contents and Storage Conditions***

Materials	Format	Quantity
Gene-specific siRNA duplexes	Purified and sequence-verified	2nmol per vial. Three unique duplexes per gene.
SR30004-Universal scrambled negative control siRNA duplex	Purified and sequence-verified	2nmol per vial.
SR30005-Rnase free siRNA duplex re-suspension buffer	100mM KAc/30 mM HEPES pH 7.5	2ml per vial

## ***Storage Conditions***

The dried duplexes can be stored at 4°C. However, once reconstituted with dH<sub>2</sub>O, the plasmids must be stored at -20°C.

## ***Related products***

- Cy3-labeled fluorescent transfection control siRNA duplex
- HPRT1 Positive Control siRNA duplex
- TrueORF™ Tagged ORF clones
- HUSH™ shRNA plasmids
- siTRAN siRNA Transfection Reagent

## ***Additional materials recommended***

- Transfection reagent: Transfection reagents must be selected and optimized based on the cell type being used.
- Cell line and cell culture supplies: user preferred
- Reagents for cell lysis: user preferred
- Reagents and supplies for immunoblots.

## ***Notice to purchaser***

This product is for research use only. Use in and/or for diagnostics and therapeutics is strictly prohibited. By opening and using the product, the purchaser agrees to the following: The plasmids may not be distributed, resold, modified for resale or used to manufacture commercial products without prior written approval from OriGene Technologies, Inc. If you do not agree to the above conditions, please return the UNOPENED product to OriGene Technologies, Inc. within ten (10) days of receipt for a full refund.

## RNAi Collection Overview

As a cellular defense mechanism, host cells process double-stranded RNA into small molecules which target homologous RNAs for destruction (Hannon 2002). In mammalian cells, RNA interference (RNAi) can be triggered by siRNAs that cause strong, yet transient inhibition of gene expression on specific genes (Elbashir 2001). These siRNAs can be synthesized and transfected into mammalian cells, resulting in effective suppression of gene expression. Unfortunately, such suppression is transient. By contrast, short hairpin RNAs (shRNA) can suppress gene expression over a prolonged period by continually expressing an RNA duplex (Brummelkamp 2002; Paddison 2002).

In cells, small interfering RNAs (siRNAs) are produced by enzymatic cleavage of long dsRNAs by the RNase-III class endoribonuclease Dicer. The siRNAs associate with the RNA Induced Silencing Complex (RISC) in a process that is facilitated by Dicer. Dicer-Substrate RNAi methods take advantage of the link between Dicer and RISC loading that occurs when RNAs are processed by Dicer. Traditional 21-mer siRNAs are chemically synthesized RNA duplexes that mimic Dicer products and bypass the need for Dicer processing. Dicer-Substrate RNAs are chemically synthesized 27-mer RNA duplexes that are optimized for Dicer processing and show increased potency when compared with 21-mer duplexes [1, 2].

The siRNA-27 kit contains three Dicer-Substrate 27-mer duplexes targeting a specific gene that are selected from a predesigned set of duplexes from the RefSeq collection of Genbank. OriGene's siRNA collection of 27-mer sites were chosen by a rational design algorithm that integrates both traditional 21-mer siRNA design rules as well as new 27-mer design criteria. In addition, analysis was performed to ensure that the chosen sites do not target alternatively spliced exons and also do not include known SNPs; these sequences are therefore optimized at several levels.

In addition to three target-specific duplexes, the siRNA kit contains a scrambled universal negative control RNA duplex (DS Scrambled Neg) that is absent in human, mouse, and rat genomes. OriGene guarantees that at two of the three Dicer-Substrate duplexes in the siRNA-27 kit will give at least 70% knockdown of the target mRNA when used at 10 nM concentration and assayed by quantitative RT-PCR when the fluorescent transfection control duplex indicates that >90% of the cells have been transfected and the HPRT positive control works with the expected efficiency.

Positive control Dicer-Substrate RNA duplex (HPRT1 Positive Control) which targets a site in the HPRT (hypoxanthine guanine phosphoribosyltransferase 1) and a Cy3<sup>TM</sup> - labeled transfection control RNA duplex are available for purchase.

## Product Application

### *Reconstitution of siRNA duplexes in duplex buffer*

1. Materials should only be handled with gloves under RNase-free conditions.
2. Briefly centrifuge each tube to ensure that all material is in the bottom of the tube and not in the cap before opening for the first time. Dried oligo often dislodges during shipping and can be lost.
3. Resuspend duplexes in the provided Duplex Buffer. Other sterile, RNase-free

buffers are suitable that contain at least 50 mM NaCl or KCl. Oligos require the presence of cations to form stable duplexes, so resuspension in water is NOT recommended. If your buffer contains Mg salts (OriGene's Duplex Buffer does not contain Mg), do not heat above 75 °C as this may result in chemical degradation of the RNA.

- a. 2 nmoles of each target specific duplex is provided. Addition of 100 µl of RNase-free Duplex Buffer will result in 20 µM final concentration; vortex thoroughly and microfuge prior to use.
- b. 1 nmole of each control duplex is provided. Addition of 50 µl of RNase-free Duplex Buffer will result in 20 µM final concentration; vortex thoroughly and microfuge prior to use.
- c. Heat to 94°C for 2 minutes, remove from heat and allow tube to cool to room temperature. The oligos were dried in duplex form so heating may not be necessary, however following this protocol ensures that the contents will be fully duplexed.
- d. Each 1 nmole tube contains enough duplex for the number of transfections shown in Table 1. Target-specific duplexes are provided with 2 nmole yield and will accordingly yield more transfections.

**Table 1. Transfections per1 nmole duplex RNA.**

Plate Format	Media Volume (µL)	# of Transfections at 1 nM	# of Transfections at 10 nM
6 well	2500	400	40
12 well	1200	820	82
24 well	600	1660	166
48 well	250	4000	400
96 well	150	6660	666

4. Once hydrated, duplexes should be stored at -20 °C or -80 °C. While generally stable to freeze/thaw cycles, OriGene recommends that daughter aliquots be made for routine use to minimize the frequency of freeze/thaw events for primary stock tubes. Minimize light exposure for dye-labeled duplexes.
5. Refills of gene-specific or control duplexes can be ordered on [www.origene.com](http://www.origene.com). Note that additional control duplexes are available. These are not included in the kit but may be useful for certain applications. Sequences can be ordered in larger scale as needed.

### **Outline of an RNAi Experiment**

1. Establish optimal transfection method for your cell type and culture media (use fluorescent-labeled transfection control duplex); favored approach is fluorescence microscopy. Greater than 90% of cells should show dye uptake when examined 4-24 hours after transfection.
2. Demonstrate that RNAi is working using positive control (HPRT1 Positive Control duplex); favored approach is quantitative real-time RT-PCR. HPRT should show >90% knockdown 24 hours post-transfection at 10 nM dose.
3. Test target specific duplexes and perform dose response curve. OriGene recommends testing duplexes at 10 nM, 1 nM, and 0.1 nM concentrations. Knockdown of mRNA levels should be assayed at 24-48 hours post transfection.

To limit off-target effects, routine studies should subsequently be performed using the lowest concentration of RNA duplex that achieves the desired level of suppression of the target mRNA.

4. Perform RNAi studies using duplexes identified as “effective by >70% reduction in RNA levels”. OriGene recommends that the results of two duplexes against the same target be compared to control for potential off-target effects and other artifacts.
  - a. mRNA levels can generally be measured 24-48 hours post transfection.
  - b. Protein levels can generally be measured at 48-72 hours post transfection, however this may vary depending on the half-life of the protein studied and cell growth rate.
  - c. Phenotype studies should parallel protein evaluation.
5. Controls: While examination of non-transfection and mock-transfection cultures (lipid or electroporation alone) are useful, OriGene recommends that control cultures transfected using control RNA duplexes be used for target level normalization. A randomized sequence (Universal Scrambled Negative Control) duplex is provided for this purpose, which is not present in human, mouse, or rat. Duplexes targeting reporter genes such as EGFP or Firefly Luciferase can also be employed as controls, however these duplexes have not been optimized for this use.

## ***Choice of Transfection Method***

A successful RNAi experiment starts with good transfection. Unfortunately, the same methods that work well for plasmid DNA transfection may not work well when using short duplex RNAs; reagents and protocols must be optimized for each type of nucleic acid transfected and for each different cell type employed. The Cy3-labeled Transfection Control and HPRT1Positive Control duplexes are available and can be purchased to assist with this phase of the RNAi project.

Cationic lipids are commonly employed to transfect duplex RNAs into tissue culture cells. It may be necessary to test a variety of lipid agents with each cell type used to identify the optimal reagent that introduces duplex RNA into cells with high efficiency yet results in minimal toxicity. After successful transfection of a fluorescent duplex, the cells will show both punctuate and diffuse cytoplasmic staining and also remain “healthy appearing” when examined using fluorescence microscopy. A Cy3-labeled transfection control duplex (excitation max 556 nm, emission max 570 nm) is available for purchase. OriGene recommends use of the fluorescent transfection control oligos at 10 nM concentration; optimal amounts may vary with the cell types and equipment employed. Transfection efficiency should be assessed within a 4-24 hour window of transfection.

Some cell types are resistant to transfection using cation lipids. Introduction of RNAi duplexes into such cells can often be achieved using electroporation. Liposomal and peptide-based transfection reagents are also available and may be of benefit. Like cationic lipid transfection, electroporation requires optimization. Important variables include voltage, pulse length, pulse number, buffer composition and RNA concentration. A recently published report suggest that the off-target effects (alteration in gene expression profiles not specifically targeted by the siRNA) can be reduced using electroporation mediated transfections as compared to those performed with lipid reagents [3].

## ***Transfection Optimization and protocol***

Although the transfection protocol below has been shown to result in highly efficient transfection, we recommend that the transfection conditions be optimized for each distinct cell line. The following variables should be considered:

- A. Cell density (% confluence at transfection): The recommended confluency for most adherent cell types at transfection is 50-70%. The optimal confluency may vary significantly among different cell lines.
- B. siTran 1.0 to siRNA ratio: We suggest that the conditions should be optimized for each cell line using the conditions in table 1 as a starting point.

Table 1. Recommended starting transfection conditions for siRNA using siTran 1.0

Plate Type	final Vol. (ul)	Vol. siRNA (5.0 uM), (10 nM final)	siTran1.0	Opti-MEM
96 well	120	0.3 ul	1.5 ul	8.5 ul
48 well	340	0.6 ul	3.0 ul	17 ul
24 well	680	1.20 ul	6.0 ul	34 ul
12 well	1200	2.0 ul	10.0 ul	60 ul
6 well	2800	4.0 ul	20.0 ul	120 ul

Table 2. 2 Recommended starting transfection conditions for siRNA and plasmid DNA in a co-transfection using siTran 1.0

Plate Type	final Vol. (ul)	Vol. siRNA (5.0 uM), (10 nM final)	Plasmid DNA	siTran1.0	Opti-MEM
96 well	120	0.3 ul	30 ng	1.5 ul	8.5 ul
48 well	340	0.6 ul	60 ng	3.0 ul	17 ul
24 well	680	1.20 ul	120 ng	6.0 ul	34 ul
12 well	1200	2.0 ul	180 ng	10.0 ul	60 ul
6 well	2800	4.0 ul	360 ng	20.0 ul	120 ul

### Protocol for transient transfection in a 96 well plate

1. Plate cells the day before the transfection. Plate the appropriate dilution of cells to provide a confluency of 50-70% in 24 hrs.
2. Dilute siRNA to 5 uM using sterilized duplex buffer supplied by siRNA manufacturer.
3. Set up a series of sterilized 0.2 ml PCR tubes. Add 10 ul of Opti-MEM and 0.3 ul of 5 uM siRNA to each well and also 50 ng plasmid DNA if performing co-transfection.

4. Prepare adequate diluted siTran 1.0 (10 ul per well) in Opti-MEM using the recommended ratio siTran (1.5ul siTran/8.5 ul Opti-MEM) and/or ratios suitable for your cell line. Add 10 ul of the siTran dilution to each PCR tube with siRNA solution.
5. Mix the solution and incubate at room temperature for 10 mins.
6. Transfer the siTran/siRNA mixtures to each well of cells in the microtiter plate
7. Incubate the plate at 37oC in a CO2 incubator.
8. View cells under a fluorescence microscope or harvest cells 24 to 48 hrs after transfection.

## Positive and Negative Controls

Fluorescent-labeled duplexes permit rapid optimization of transfection protocols. The best control for a good transfection, however, is demonstration of successful suppression of a control gene using a known effective siRNA. OriGene offers a validated HPRT1 Positive control duplex for this purpose. It is a potent Dicer-Substrate duplex that has been shown to suppress HPRT mRNA levels >90% when used at 10 nM concentration. OriGene also provides a validated negative control duplex as part of the kit that is not present in the human, mouse, or rat genomes.

Validated primers for use in a HPRT SYBR Green Q-PCR assay are available from OriGene's qSTAR primer collection.

## Measuring RNAi Knockdown of Gene Expression

OriGene recommends using quantitative assessment of mRNA levels as the primary tool to establish efficacy of RNAi methods. Measurements of protein levels, enzyme activity, or phenotype can vary widely depending on protein half-life, turnover rates, and other factors such that it is possible to observe a seemingly negative result even when mRNA levels have been substantially suppressed. Protein levels and phenotype effects can be assessed after it has been established that mRNA levels are reduced.

**RNA** □□□  $\implies$  **Protein**  $\implies$  □□□**Phenotype**

## FAQ: siRNA-27 kit

### What are Dicer-Substrate RNAs?

Dicer-Substrate RNAs are chemically synthesized RNA duplexes that are optimized for Dicer processing and show increased potency when compared with 21-mer duplexes.

### What is the difference between DsiRNAs (27 mers) and traditional 21mer siRNAs?

Many researchers today employ synthetic 21-mer RNA duplexes as their RNAi reagents, which mimic the natural siRNAs that are the result of Dicer processing of long substrate RNAs. An alternative approach is to use synthetic RNA duplexes (DsiRNAs) that are greater than 21-mer in length, which are substrates for Dicer. DsiRNAs are processed by Dicer into 21-mer siRNAs and designed so that cleavage results in a single, desired product. This DsiRNA design provides Dicer with a single favorable PAZ binding site that helps direct the cleavage reaction. Functional polarity is introduced by this processing event, which favors AS strand loading into RISC, and the increased potency of these reagents is thought to relate to increased DsiRNAs being loaded into and processed

through Dicer. The Dicer-substrate approach can result in reagents having as much as 10-fold higher potency than traditional 21-mer siRNAs at the same site.

### **Is the designed RNA duplex specific?**

Our Bioinformatics department did a full screening of the recommended sequences before adding them to the database. This included blasting the sequences to ensure they only line up with the gene in questions, performing secondary structure checks, and making sure the end stabilities were favorable for incorporation into the RISC complex.

### **How do I order GFP controls to go along with the siRNA kit?**

We do not offer any GFP controls along with the siRNA kit. The controls we do offer are listed below

Cy3-labeled transfection control RNA duplex (SR30002)

HPRT1 Positive Control duplex (SR30003)

Universal scrambled negative control duplex (SR30004)

These duplexes can be ordered separately at [www.origene.com/siRNAcontrols](http://www.origene.com/siRNAcontrols)

### **Does the siRNA kit come with primers and RT-PCR materials?**

The siRNA kit includes:

3 Dicer-Substrate RNAi duplexes, 2nmoles each- to test knockdown you will need to order primers separately.

1 negative control duplex

Rnase-free duplex re-suspension buffer

qSTAR primer pairs and SYBR Green master mix are available for purchase. For more information, visit [www.origene.com/geneexpression/](http://www.origene.com/geneexpression/)

### **What concentration of DsiRNAs do I need to use in my experiment to see knockdown?**

For an RNAi or antisense experiment, the actual level of target gene knockdown relates to the transfection efficiency. A positive control such as HPRT siRNA should always be used in each experiment to assess transfection efficiency.

In addition, varying amounts of test DsiRNA ranging from 0.1 to 50 nM can be used to determine which DsiRNA shows maximum knockdown

### **Do you provide any non-silencing controls for Dicer Substrate RNAi duplexes?**

We currently include a scrambled negative control for our siRNA kit. Other controls are available for purchase separately.

The normal turn around time for the siRNA kit is 7-9 business days. The number of assays that can be done will vary widely with the transfection reagent used and the cell line you are trying to transfect. Optimization will be necessary. Normal optimization concentrations vary between 1-50nM siRNA, with 10nM being most common.

### **What is OriGene's guarantee on siRNA duplexes?**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will give at least 70% or more knockdown of the target mRNA when used at 10 nM concentration and assayed by quantitative RT-PCR when the fluorescent transfection control duplex indicates that >90% of the cells have been transfected and the HPRT positive control works with the expected efficiency.

## References

1. Kim, D.H., et al., Synthetic dsRNA Dicer substrates enhance RNAi potency and efficacy. *Nat Biotechnol*, 2005. **23**(2): p. 222-6.
2. Rose, S.D., et al., Functional polarity is introduced by Dicer processing of short substrate RNAs. *Nucleic Acids Res*, 2005. **33**(13): p. 4140-56.
3. Fedorov, Y., et al., Different delivery methods-different expression profiles. *Nat Methods*, 2005. **2**(4): p. 241.
4. Amarzguioui, M., Improved siRNA-mediated silencing in refractory adherent cell lines by detachment and transfection in suspension. *Biotechniques*, 2004. **36**(5): p. 766-8, 770.
5. Bieche, I., et al., Quantitation of hTERT gene expression in sporadic breast tumors with a real-time reverse transcription-polymerase chain reaction assay. *Clin Cancer Res*, 2000. **6**(2): p. 452-9.



### UK & Rest of World

184 Millon Park, Abingdon  
OX14 4SE, Oxon, UK  
Tel: +44 (0) 1235 828 200  
Fax: +44 (0) 1235 820 482

### Switzerland

Centra Nord-Sud 2E  
CH-6934 Bioggio-Lugano  
Tel: +41 (0) 91 604 55 22  
Fax: +41 (0) 91 605 17 85

### Deutschland

Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
Tel: +49 (0) 69 7790999  
Fax: +49 (0) 69 13376880

### North America

23591 El Toro Rd, Suite #180  
Lake Forest, CA 92630  
Tel: +1 800 987 0985  
Fax: +1 949 265 7703

amsbio

info@amsbio.com

www.amsbio.com  
AMS Biotechnology