




Oscilloscope Dual Trace/Channel

241-541



Instruction Manual

Catalogue

1. Brief Introduction

2. Technical Specifications

3. Matters needing paying attention to

3.1 Opening

3.2 Voltage & Power Checking

3.3 Environment

3.4 Installation & Operation

3.5 Fluorescent Coating for CRT

3.6 Max. Voltage on Input Terminal

4. Operation Methods

4.1 Description on Front Panel

4.2 Description on Rear Panel

4.3 Basic operation & Single-channel Operation

4.4 Dual-channel Operation

4.5 Plus-minus Operation

4.6 Selection of Trigger Source

4.7 Sweep-speed controlling

4.8 Sweep Magnification

4.9 X-Y Operation

4.10 Probe Calibration

4.11 Component Test (only for specifically model)

5. Measurement

5.1 Checking & Adjusting before measuring

5.2 Amplitude Measuring

5.3 Time Measuring

5.4 TV Signal Measuring

6. Maintenance

6.1 Fuse Replacement

6.2 Cleaning

6.3 Standard layout

1. Brief Introduction

Characterized in beautiful appearance, reasonable internal structure, advanced functions and lower price, the serial dual-trace oscilloscope has been widely used in the fields of teaching, enterprises, scientific research and medical treatments. The main functions and characteristics are listed as following:

- The vertical frequency response is up to 20MHz / 30MHz and the maximum sensitivity would be up to 1mV/DIV.
- The sensitivity would be 20V/DIV with the broad voltage input range.
- The vertical operating modes would be CH1, CH2, ALT, CHOP, CH1+CH2, CH1-CH2.
- There are three trigger sources for selection: INT, EXT, Power (ALT function for INT trigger mode).
- The triggering and sweeping modes are: NORM, AUTO (TV-H) , TV-V & Level Lock.
- The lowest sweep time would be 0.5s/DIV, the fastest would be 0.2 us/DIV. The sweep speed can be up to 20ns/DIV by ×10 magnification.
- The sweep can be adjusted with hold-off time and it can synchronize the periodic and complex waveforms.
- It has the function of X-Y display and the highest sensitivity would be 1mV/DIV.
- The auto level-locking function can display waveforms stably without adjusting the level potential.
- The light emitting diode (LED) would be lighted on when the stable waveforms are displayed.
- The auto-focusing function can display waveforms clearly when the CRT lightness is changed.
- It has the function for components test (only for specifically model).
- It has DDS function signal output (only for specifically model).

2. Technical Specifications

Specification	Items	20MHz Oscilloscope	30MHz Oscilloscope
---------------	-------	--------------------	--------------------

		20MHz Oscilloscope	30MHz Oscilloscope
Vertical System	Sensitivity	5mV~20V/DIV in 1-2-5 sequence, altogether 12 steps, CH1, CH2 to 1mV/DIV with ×5 magnification	
	Accuracy	×1: ≤±3%, ×5 MAG: ≤±5%	
	Variable Ratio	≥2.5:1	
	Bandwidth (-3dB)	×1: DC(AC10Hz)~20MHz, ×5: DC(AC10Hz)~7MHz	×1: DC(AC10Hz)~30MHz, ×5: DC(AC10Hz)~7MHz
	Rising Time	×1: ≤17.5ns, ×5: ≤50ns	×1: ≤12ns, ×5: ≤50ns

	Input Impedance	1MΩ±5% // 25pF±5pF		
	DC Balance	5mV~20V/DIV: ±0.5DIV		
	Linearity	The amplitude change would be within±0.1V when the waveform moves vertically in the middle of the division.		
	Vertical Mode	CH1, CH2, ALT, CHOP, ADD (CH1+CH2, CH1-CH2)		
	Input Coupling	AC, GND, DC		
	Max. Input Voltage	400V with the frequency ≤1kHz Max. effective readout would be 160Vp-p (56Vrms sine wave) when the probe is set as 1:1. Max. effective readout would be 400Vp-p (140Vrms sine wave) when the probe is set as 10:1.		
	CH2 INV BAL	≤1DIV		
Trigger	Trigger Sources	INT, EXT, LINE INT		
	Trigger Source	CH1, CH2, VERT.		
	Trigger Modes	NORM, AUTO (TV-H), TV-V, LEVEL LOCK		
	Coupling	AC: 5Hz to the whole frequency range		
	Polarity	+/-		
	Sensitivity	INT: 5Hz~10MHz≤1DIV; 10MHz~20MHz≤1.5DIV; ≤2DIV	INT: 5Hz~10MHz≤1DIV; 10MHz~30MHz≤2DIV; ≤2DIV	
		EXT: 5Hz~10MHz≤200mVp-p; 10MHz~20MHz≤300mVp-p; ≤500mVp-p	EXT: 5Hz~10MHz≤200mVp-p; 10MHz~30MHz≤400mVp-p; TV: ≤500mVp-p	
Input impedance with EXT trigger signals Max Input Voltage	1MΩ±5%//25pF±5pF 400V(DC+ACpeak) AC frequency:≤1kHz			
Sweep time	0.5s~0.2us/DIV, in 1-2-5 sequence, altogether 20 steps			
Accuracy	×1: ≤±3%; ×10MAG: ≤±5% (20ns~50ns : ±10%)			
Variable Ratio	≥2.5:1			
Linearity	×1: 5%; ×10MAG: 10% (20ns~50ns : 15%)			
Movement by ×10 MAG	<2DIV in CRT center			

Specification		Items	20MHz Oscilloscope	30MHz Oscilloscope
		X-Y Mode	Sensitivity	Same the vertical system
Frequency Bandwidth (-3dB)	DC: 0~500kHz; AC: 10Hz~500kHz			
X-Y Phase Difference	≤3°(DC~50kHz)			
Calib rating Signa	Waveform	Square wave		
	Frequency	Approx. 1kHz		

	Output voltage	2Vp-p±2%
	Output Resistance	Approx. 1kΩ
CRT	Model	15SJ118Y14
	Color & Afterglow	Green, middle
	Effective Screen Area	8×10DIV [1DIV=10mm(0.39in)]
	Scale	Internal
	Trace Rotation	Adjustable on panel

Power Requirement:
Voltage: Fixed AC220V±10%
50Hz/60Hz
Power Consumption: about 40VA

Operation Environment:
Used indoors Frequency:
Sea level: 2000m
Environment temperature: 0~40℃
Humidity: 85% RH, Dry
Size: 310W×145H×440D(mm)
Weight: Approx. 8kg
Storage Temperature:-10~70℃

Z Axis (Only for specifically model, other models optional)

Item	Specification
Sensitivity	5V p-p (Trace becomes brighter with negative input)
Bandwidth	DC-2MHz Input
Impedance	Approx. 10kΩ
Maximum Input Voltage	30V (DC+AC peak) AC frequency:≤1kHz

CH1 Signal Output (Only for specifically model, other models optional)

Item	Specification
Output Voltage	Least 20mV/DIV
Output Impedance	Approx. 50Ω
Bandwidth	50Hz-5MHz (-3dB)

DDS function signal output (Only for specifically model)

Item	Specification
Frequency Range	0.2Hz~2MHz
WaveForm	Sine,Square, Triangle Wave
Frequency Variable Range	10:1 or more
Voltage	Over 20 V p-p

DC offset	$\pm 10V$ when open
Sinewave Distortion	0.3%
Squarewave Unsynnetrr	$\pm 3\%$ or less
Squarewave Rise/Fall Time	50 Ω output:30ns or less

3. Matters needing paying attention to

3.1 Opening

The oscilloscope has been inspected and tested strictly before it is taken out of the factory. Please check whether it is damaged during delivery after receiving it. And please contact with the supplier or the consigner if there is any problem with the machine.

3.2 Voltage & Power Checking

Please check whether the voltage used can meet the voltage requirement before the machine is put to use.

Pay attention: The wrong voltage used would damage the oscilloscope.

!!! Warning: The ground protection terminal should be well connected to avoid the electric shock.

Please replace the specific fuse listed in the following table when it is burnt.

Power voltage	Range	Fuse
AC220V	198~242V	1A

!!! Warning: The power should be shut off while replacing fuses to protect human beings.

3.3 Environment

The environment temperature would be in the range of 0~40 $^{\circ}C$. Otherwise the circuits inside the machine would be damaged.

3.4 Installation & Operation

Please be sure that the holes for heat radiating on the oscilloscope would not be plugged.

The auto-protection function would be weakened if the machine is not used under the specific condition.

3.5 Fluorescent Coating for CRT

Please don't set the CRT brightness to be the brightest state or leave the dot to be in one position for a long time in order that the CRT fluorescent coating can be protected.

3.6 Max. Voltage on Input Terminal

The following table lists the maximum voltage on the input terminal and the probe. When the probe is set as 1:1, the maximum effective readout voltage is 160Vp-p(56Vrms at the sine wave). If the probe is set as 10:1, the maximum effective readout voltage is 400Vp-p(140Vrms at the sine wave).

Input Terminal	Max. Input Voltage
----------------	--------------------

CH1, CH2	400Vp-p
Ext Trigger In	400Vp-p
Probe	400Vp-p

!!! Warning: Don't use it excess the above values to protect the machine. The frequency of the maximum input voltage would be less than 1kHz.

Fig. 4-1A




Fig. 4-1B




Fig. 4-2A





Fig. 4-2B



4. Operation Methods

4.1 Description on Front Panel (Refer to Fig. 4-1A or Fig. 4-1B)

CRT:

- ⑦ —Power: Main power switch. The LED ⑥ would be lighted When it is on.
- ① —Intensity: Adjust the intensity of the trace or the dot.
- ③ —Focus: Adjust the clearness of the trace or the dot.
- ④ —Trace Rotation: The half-fixed potentiometer is used to adjust the horizontal trade to be parallel with the scale.
- ③⑩ — Color Filter: Make the display effect be more comfortable.

Vertical Axis:

- ①⑦ —CH1 (X) Input: Channel 1 input terminal, which is used as X axis input in X-Y mode.
- ①⑧ —CH2 (Y) Input: Channel 2 input terminal, which is used as Y axis input in X-Y mode.
- ②⑧ ③③ — ×5 MAG for CH1 & CH2: Used to adjust the vertical sensitivity to be 1mV/DIV.
- 15 16 —AC-GND-DC: Select the input mode for the input signals of the vertical axis.



AC: AC coupling

GND: The input of the vertical magnifier is grounded and the input is cut off.

DC: DC coupling

⑪ ⑫ —Vertical Attenuator: Adjust the vertical deflection sensitivity from 5mV to 20V/DIV in 12 steps.

⑬ ⑭ —Vertical Variable: The variable ratio would be no less than 2.5:1. The sensitivity calibration would be the label value in the calibrated position.

⑧ ⑨ --▲▼ Vertical Position: Adjust the vertical position of the trace on the screen.

⑩ —Vertical Mode: Select the operating modes for CH1 & CH2 magnifiers.

CH1 or CH2: Channel 1 or 2 displayed separately.

Dual: Two channels are displayed at the same time.

⑳ —Dual-trace Display:

ADD: Press to display the algebraic sum of the two channels CH1+CH2 in the dual-trace display mode. Press CH2 INV 35 , the display would be the algebraic difference CH1-CH2.

ALT: Press to display the traces alternatively of CH1 and CH2 in the dual-trace display mode (usually for the quick-speed sweep).

CHOP: Press to display the traces on CH1 & CH2 in Chop mode.

㉑ —CH2 INV: The signal in CH2 would be inversed. Press the key to make the signals from CH2 and the inner trigger signals of CH2 be inversed at the same time.

Trigger:

㉒ —EXT trigger input: For the external trigger signals. The trigger source selector should be set to be EXT when it functions. (Only for specifically model)

㉓ —Trigger Source Selection:

INT: Select the signals from CH1 or CH2 as the trigger source.

LINE: Select AC power as the trigger signal.

EXT: The external signals on ㉒ is selected as the trigger signal.

㉔ —INT Trigger source selection:

VERT: When the vertical mode selector ⑩ is set to be DUAL and the trigger source selector ㉓ is set to be INT, press VERT to select CH1 & CH2 alternatively as the internal trigger signal source.

CH1: Select Channel 1 as the internal trigger signal source.

CH2: Select Channel 2 as the internal trigger signal source.

㉕ —Polarity: Select the polarity of the trigger signals. "+" means the rising-edge trigger and "-" means the trailing-edge trigger.

②④—Trigger Level: Display a synchronized stable waveform and set the starting point for the waveform. Rotate clockwise to increase the level and rotate counterclockwise to decrease the level.

②⑦—Trigger mode: Select the trigger mode.

AUTO (TV-H): The sweep mode is the AUTO mode when there is no trigger signal input.

Used to observe the TV-H signals. (It can only be synchronized when the synchronized signal is of negative pulse.)

NORM: The trace wouldn't be displayed when there is no trigger signal.

TV-V: Used to observe the TV-V signals.

(It can only be synchronized when the synchronized signal is of negative pulse.)

LOCK: The trigger lock is locked in a fixed level. It is not necessary to adjust the level to get the synchronized signal when the sweep speed or the signal amplitude is

changed. ③④ —Trigger Indicator: The LED would be lighted in the trigger sweep.

Time base:

②⑥—Horizontal Sweep Speed Switch: The sweep speed would be divided into 20 steps from 0.2us/DIV to 0.5s/DIV. (Only for specifically model)

④③—Horizontal Sweep Speed Switch. (Only for specifically model)

The sweep speed would be divided into 20 steps from 0.2us/DIV to 0.5s/DIV.

X-Y: X-Y mode.

②⑤—Horizontal Variable: Adjust the horizontal sweep time to calibrate it to be the same shown by TIME/DIV on the panel. The TIME/DIV can be adjusted continuously and would be in calibrated position when it is rotated to the end clockwise. The whole time delay can be up to 2.5 times or more.

② —HOLD OFF: Change the sweep off time, synchronize multicycle waves. (Only for specifically model)

②⑨- ◀▶ Horizontal Position: Adjust the horizontal position of the traces on the screen.

③①—Sweep Mode: Select the sweep mode. (Only for specifically model)

×1: The sweep is not magnified.

×10: The sweep is 10 times magnified.

X-Y: The key ②⑥ doesn't function when in X-Y mode.

④②—Sweep Mode: Select the sweep mode. (Only for specifically model)

×1: The sweep is not magnified.

×10: The sweep is 10 times magnified.

DDS function:

- ④④--Voltage and DC offset regulation button
- ④⑤--DC offset:when open DC regulation
- ④⑥--Adjust the frequency knob
- ④⑦--Stall frequency adjustment button
- ④⑧--The frequency of stalls instructions
- ④⑨--Waveform select button
- ⑤①--Waveform selection of indicators
- ⑤①--Signal output

Other keys:

- ⑤-- CAL: Supply a square wave signal with the amplitude 2Vp-p and the frequency 1kHz to calibrate the compensate capacity of the 10:1 probe and detect the deflection factors of the horizontal and vertical systems of the oscilloscope.
- ①⑨--GND: The ground terminal for the oscilloscope.
- ③⑥--Input Hole for Component testing(only for specifically model).
- ③⑦--Oscilloscope/Component test: Press it to test components. (only for specifically model).

4.2 Description on Rear Panel (Refer to Fig. 4-2A or Fig. 4-2B)

- ③⑧--Power plug & Fuse: 220V voltage plug.
- ③⑨--Z-Axis Input: Input terminal for external intensity modulation signal. (option)
- ④①--Signal Output: CH1 signal output terminal,Used to measure frequency. (option)
- ④①--EXT trigger input: For the external trigger signals. The trigger source selector should be set to be EXT when it functions. (Only for specifically model)

4.3 Basic operation & Single-channel Operation

Check whether the voltage connected meets the requirements of the machine and set the following keys listed in the table below:

Function	No.	Setting
POWER	7	OFF INTEN
	1	MIDDLE
FOCUS	3	MIDDLE
VERT MODE	10	CH1 DUAL-
TRACE DISPLAY	32	ALT
CH2 INV	35	Release
VERT POSITION	8,9	MIDDLE
VOLTS/DIV	11,12	50mV/DIV

VAR	13,14	CAL
AC-GND-DC	15,16	GND
TRIGGER SOURCE	21	INT
SLOPE	23	+
INT TRIGGER SELECTOR	22	CH1
TRIGGER MODE	27	AUTO
TIME/DIV	26	0.5ms/DIV
SWP.VAR	25	CAL HOR.
POSITION	29	MIDDLE
SWEEP MODE	31	×1

After finishing the above setting, connect the power cable and:

- (1) Turn on the power. The trace would appear on the screen about 20 seconds later. If the trace wouldn't appear after 60 seconds, please check the setting of the switches and controlling knobs again.
- (2) Adjust the knobs INTEN and FOCUS to get the clear and comfortable trace.
- (3) Adjust CH1 position and the potentiometer for trace rotating to make the trace be parallel with the horizontal scale. (Use a screw to rotate the potentiometer ④).
- (4) Use the probe 10:1 to input the calibrating signals into CH1 input.
- (5) Set the key AC-GND-DC to be AC. A square wave shown as Fig. 4-3 would appear on the screen.
- (6) Adjust the FOCUS to make the trace be clear.
- (7) As for other kinds of waveforms, adjust the vertical attenuator, the sweep time switch, the vertical and horizontal position knob to the proper setting. The amplitude and the time can be read out easily.

The operation of CH2 is the same as that of CH1.




Fig4-3




Fig4-4

4.4 Dual-channel Operation

Set the vertical mode to be DUAL mode and DUAL display mode is set to be ALT. the trace on CH2 would display on the screen. CH1 would display a square wave (from the calibrating signal) and CH2 would display a straight line (because no signals are input to the channel). Connect the calibrating signal to CH2 input and set AC-GND-DC to be AC mode. Adjust the vertical position Ⓢ & Ⓢ to get the waveforms shown as Fig. 4-4. The signals from CH1 & CH2 would be displayed alternatively on the screen, which can be used to observe the signals from two ways with short time sweeping. When the DUAL display is setting to be CHOP, the signals from CH1 & CH2 would be displayed on the screen separately with the speed of 250kHz, which can be used to observe the signals from two ways with long time sweeping. In two channels' operation, select the signals from CH1 or CH2 as the trigger signals by the trigger source selector in the mode of DUAL. If the signals from CH1 &* CH2 are synchronized, the two waveforms can be displayed stably. Or one channel is for one stable wave. If the internal trigger source selector 22 is set as VERT, the two waves can be stably displayed.

4.5 Plus-minus Operation

Set DUAL trace display 32 to be ADD by setting Vertical mode 10 to be DUAL and the algebraic sum of CH1 & CH2 can be displayed. If CH2 INV is pressed, the algebraic difference would be displayed. And at this moment, the attenuators of the two channels should be set the same. The vertical position can be adjusted by the knob ▲ ▼ Vertical Position. Set the knob to be in the middle since there is a linear change on the vertical magnifier.

4.6 Selection of Trigger Source

It is very important for you to select the effective trigger sources. You should be very

familiar with the selection, the functions and the operation of the trigger source.

(1) Trigger mode selector:

AUTO (TV-H): The sweep generator would produce a sweep signal by free vibration. If there is a trigger signal, it would change to the trigger mode automatically. Set it to be AUTO when the first wave is observed. Set it to the proper position after a stable wave is observed. Reset the mode to be AUTO after Setting other knobs as specified. If DC signal or small signal is measured, AUTO mode must be used. Set the mode to be AUTO (TV-H) to observe TV-H signals. The sweep time should be set to be $10 \mu\text{s}/\text{DIV}$. Demonstrated that several line of waveforms, may use the vernier adjustment knob adjustment trace time to arrive need the number of lines. Sends in oscilloscope's synchronized signal to be the negative polarity. shown as Fig. 4-5:




Fig. 4-5

NORM: There is no display on the screen when the sweep generator is in the statistic mode. Sweep for one time when the trigger signal passes through the gate level set by the trigger level switch, then the sweep generator returns back the statistic state. One trigger is for one sweep. When the dual-trace display is set to be ALT or CHOP, there is no display except there is sufficient trigger level from CH1 and CH2.

TV-V: Set the mode to be TV-V to observe TV-V signals. The sweep time should be set to be $2\text{ms}/\text{DIV}$ (one-frame signal) or $5\text{ms}/\text{DIV}$ (interlacing sweep signal of two frames) to synchronize TV signals.

In order to get a stable wave, one signal related to the display signal in time should supplied to the trigger circuit. The trigger source selector is used to select the trigger signals.

CH1: Internal trigger mode in most cases.

CH2: The signals to the vertical input would be divided and one part would be sent to the trigger circuit before it is preset. A stable wave would be displayed on the screen since the trigger signal is the signal needing measurement.

In DUAL mode, the trigger signals are selected by the trigger source selector.

LINE: Use the frequency of the AC power from the electrical net as the trigger signal. It is very effective to measure the signal relative with the power frequency, such as: AC noise of the sound system, SCR circuits and etc.

EXT: Use the external signals to generate the sweep trigger circuit. The signal should have some relation with the measured signal in time. The wave would be triggered and displayed by the external signals.

(3) Trigger Level and Slope switch

The trigger signal would pass through a gate level when it comes into being. Adjust the trigger level knob to change the level. The level would increase in + direction and decrease in – direction. And it would be the average value when it is in the middle.

The trigger level can be used to set the starting point of the wave. For the sine wave, the starting phase can be changed. Please note that if the trigger level is adjusted to be over negative or positive, there would be no sweep signals since the trigger level is excess the amplitude of the synchronized signals.

Set the slope to be + to get the rising-edge trigger. Set the slope to be – to get the trailing-edge trigger (shown as Fig. 4-6).

Trigger Level Lock:

When the trigger mode 27 is set to be LOCK, the trigger level is locked to be a fixed value. It is not necessary to adjust the trigger level to get a stable wave even if the amplitude and the frequency of the signal are changed at this moment.




Fig4-6

The function is effective when the amplitude of the input signal or the external trigger signal is in the following range:

50Hz~5MHz \geq 1DIV(EXT:0.5V)

5MHz~20MHz \geq 2DIV(EXT:1V)

(4) VERT switch

When the vertical mode is set to be DUAL ALT display, the switch is used for alternative trigger and display. One alternative trigger signal is for one sweep cycle in ALT mode, which is helpful to test the amplitude and period of the wave, also two unrelatable signals in frequency can be observed. But it is unnecessary to measure the phase and time. One synchronized signal should be used for trigger for two channels.


4.7 Sweep-speed controlling

Adjust the sweep speed knob to select the number of the waves you want to observe. If there are too many waves on the screen, you can set the sweep speed much faster.

And if there is only one wave on the screen, you can set the sweep speed much slower. If the speed is too fast, you can only observe one part of the cycle signal. If the measured signal is of square wave, the displayed on the screen would be only a straight line.

4.8 Sweep Magnification

Very fast sweep speed is needed to observe one part of a wave. If the part is far from the starting point, it maybe out of the screen. In such case, the sweep magnification switch is useful. Press the switch, the range displayed is 10 times of the original. The speed would be 1/10 of the original. For example: 1 μ s/DIV to 0.1 μ s/DIV.



Adjusts ◀▶ position knob to be possible to observe the entire region the profile

Fig.4-7

4.9 X-Y Operation

Set the sweep mode switch to be X-Y.

X-axis: CH1 input

Y-axis: CH2 input

Note: Please pay attention to the difference of the frequency and phase between X-axis and Y-axis when the HF signal is in X-Y mode.

The oscilloscope can do so many testing in X-Y mode compared with that it works in normal mode. CRT can display an electrical picture or two transient levels which would be the direct comparison of the two levels, just like the video color-bar picture shown by the vector oscilloscope. If such kind of dynamic factors are switched into the voltage signals, their pictures can be easily displayed in X-Y mode, for example: the amplitude and the frequency, which Y-axis corresponds to the signal amplitude and X-axis corresponds to the signal frequency (Shown as Fig. 4-8).

In some cases, Lissajous figures can be observed in X-Y mode. Input the sine wave signals from X-Y input, the Lissajous figures would be displayed on the screen. The relations of the frequency and the phase between the two signals can be calculated accordingly (shown as Fig. 4-9).




Fig.4-8

Phase Difference	Display Wave			
0°				
45°				
90°				
f(y) : f(x)	1:1	2:1	3:1	3:2

Fig.4-9

4.10 Probe Calibration

The oscilloscope probe can be used for a very wide frequency range mentioned above, but the phase should be compensated. The distorted wave would lead to the error measurement. So, the probe should be calibrated before measurement (the method is referred to Paragraph 1.2 of Chapter 5).

4.11 Component Test (only for specifically model)

Set the sweep mode 31 to be X-Y mode, the input coupling switch 15 & 16 on CH1 & CH2 to be GND, CH1 attenuator 11 to be 5V/DIV, CH2 attenuator 12 to be 2V/DIV. Press the knob Oscilloscope/Component test. Insert the measured component into the hole 35. The testing clamp can also be used. The typical pictures measured are shown as Fig. 4-10:




RESISTANCE	LARGE  (20k Ω)	MIDDLE  (2k Ω)	SMALL  (200 Ω)
CAPACITANCE	LARGE  (10 μ F)	MIDDLE  (1 μ F)	SMALL  (0.01 μ F)
INDUCTANCE	LARGE  (10mH)	MIDDLE  (5mH)	SMALL  (1mH)
MANOSTAT		DICDE 	
MANOSTAT AND CAPACITANCE PARALLEL CONNECTION			RESISTANCE AND MANOSTAT IN SERIES 

Fig. 4-10

5. Measurement

5.1 Checking & Adjusting before measuring

The following items should be rechecked to keep the correct measurement and high accuracy before measurement.

5.1.1 Trace Rotation

The horizontal trace displayed on the screen would be parallel with the horizontal scale in normal cases. But there would be a slight incline on the horizontal trace because of the earth magnetic field or some other factors. So, you should check and examine the machine as following before using:

- (1) Preset the knobs on the panel to get a horizontal sweep line.
- (2) Adjust the vertical position to keep the sweep baseline on the horizontal scale on the vertical center.
- (3) Check whether the sweep baseline is parallel with the horizontal scale. If not, adjust the "Rotation" potentiometer on the front panel with a screw.

5.1.2 Probe Compensation

The probe compensation is used to compensate the error resulted from the feature difference input from the oscilloscope. The detailed procedures are listed as following:

- (1) Set the knobs on the panel (shown as Table 3) to get a sweep baseline.

- (2) Set V/DIV to be 50ms/DIV.
- (3) Connect 10:1 probe to CH1 and connect to the calibrating signal 5.
- (4) Operate the controlling knobs mentioned as Chapter 4 to get a wave like Fig. 5-1.
- (5) Observe whether the compensation is proper. If not, adjust the probe compensation component shown as Fig. 5-2.
- (6) Set the vertical mode to CH2, and connect 10:1 probe to CH2. Check and adjust CH2 probe according to the step 2 to 5.




Fig.5-1




Fig.5-2

5.2 Amplitude Measuring

5.2.1 Vp-p Measurement

Please measure the Vp-p value of the measured signal according to the following steps:

- (1) Input the signal to CH1 or CH2 and set the selected channel to VERT.
- (2) Set the voltage attenuator and observe the waveform. Keep the displayed waveform to be about 5 DIV. Check the variable and rotate to CAL clockwise.
- (3) Adjust the level to keep the waveform stable (not necessary if the level is locked)
- (4) Adjust the sweep speed to keep at least one wave cycle to be displayed on the screen.
- (5) Adjust the vertical position to keep the wave bottom to be on a certain horizontal coordinate on the screen (shown as Point A in Fig. 5-3)
- (6) Adjust the horizontal position to keep the wave top to be on a certain horizontal

coordinate in the center of the screen (shown as Point B in Fig. 5-3)

(7) Read out the division numbers between Point A & B in vertical direction.

(8) Work out the V_{p-p} of the measured signal using the following formula:

$$V_{p-p} = \text{DIV (vertical direction)} \times \text{Vertical Deflection Factor}$$

For example: in Fig. 5-3, the vertical division number between Point A & B is 4.2 DIV, the vertical deflection factor of 10:1 probe is 2V/DIV, so, $V_{p-p} = 2 \times 4.2 = 8.4V$




Fig. 5-3




Fig. 5-4

5.2.2 DC Voltage Measurement

Measure DC Voltage following the steps below:

(1) Set the knobs on the panel to get a sweep baseline.

(2) Set the coupling mode of the selected channel to be GND.

(3) Adjust the vertical position to keep the sweep baseline to be on a certain horizontal coordinate (shown as Fig. 5-4 Before Measurement) and define the point as the voltage zero value.

(4) Feed in the measured voltage to the plug of the selected channel.

(5) Set the input coupling to be DC. Adjust the voltage attenuator to keep the sweep baseline to be on a proper position on the screen. Rotate the variable clockwise to the calibrated position first.

(6) Read out the deflected division number in vertical direction, shown in Fig. 5-4, after measurement.

(7) Work out DC voltage using the following formula:

$$V = \text{DIV (vertical direction)} \times \text{Vertical Deflection Factor} \times \text{Deflection Direction}(+/-)$$

For example, in Fig. 5-4, the sweep baseline is 4 divisions above the original baseline,

the vertical deflection factor is 2V/DIV, so: $V=2 \times 4 \times (+)=+8V$

5.2.3 Amplitude Comparison

Use the following steps to measure the amplitude differences between two signals:

- (1) Feed in the reference signals into CH1 or CH2. Set the selected channel to VERT.
 - (2) Adjust the voltage attenuator and the variable to make the displayed amplitude be 5 divisions vertically, from 0% to 100% in scale.
 - (3) Keep the original setting of the voltage attenuator and the variable and switch the probe from the reference signal to the compared signals. Adjust the vertical position to make the wave bottom directly to the 0% scale on the screen.
 - (4) Adjust the horizontal position to keep the wave top on the vertical scale in the center of the screen.
 - (5) Read out the percentage (1 division equals to 4%) from the vertical coordinate in the center of the screen according to the 0% & 100% percentage standard on the left screen.
- For example: In Fig. 5-5, the dashed line indicates the reference signal and the real line indicates the compared signals. The vertical amplitude is 2DIV, so the amplitude would be 40% of the reference signal.




Fig.5-5

5.2.4 Algebra superposition

Use the following steps to measure the algebra sum or difference of two signals:

- (1) Set the vertical mode to be DUAL and select ALT or CHOP if necessary.
- (2) Feed in the two signals into CH1 & CH2.
- (3) Adjust the voltage attenuator to make the displayed amplitudes of the two signals be moderate and VOLS/DIV be the same. Adjust the vertical position to make the

waves of the two signals be in the center of the screen.

(4) Set DUAL to be ADD to get the algebra sum displayed. If you want to observe the algebra difference, please press in the knob CH2 INV. Fig. 5-6 shows the algebra sum and difference of two signals.




Fig.5-6

5.3 Time Measuring

5.3.1 Measurement of Time Intervals

The time intervals can be measured as following steps:

- (1) Feed in the signal to CH1 or CH2 and set the selected channel to be VERT.
- (2) Adjust the level to keep the wave to be displayed stably (It is unnecessary to adjust the level if it is locked).
- (3) Rotate the variable clockwise to the calibrated position. Adjust the sweep speed switch to display one to two signal cycles on the screen.
- (4) Adjust the vertical and horizontal positions to make the two points on the wave measured be on the horizontal scale in the center of the screen.
- (5) Measure the horizontal scale between the two points and work out the time interval using the following formula:

Time interval (T)=[Horizontal distance between two points (DIV)×Sweep time factor (Time/DIV)]/Horizontal Magnification Rate

For example: In Fig. 5-7, the horizontal distance between A and B is 8DIV, the sweep time factor is 2 μ s/DIV, the horizontal magnification rate is $\times 1$, so: Time Interval = (2 μ s/DIV \times 8DIV)/1=16 μ s.




Fig.5-7

5.3.2 Measurement of Cycle and Frequency

In Fig. 5-7, the measured time interval is the cycle of the signal T, the frequency would be $1/T$. For example: $T=16\mu s$, then the frequency would be : $f=1/T=1/16\times 10^{-6}=62.5\text{kHz}$

5.3.3 Measurement of Rising Time and Trailing Time

The measuring method of the rising time and trailing time is just the same as that of the time interval, only except the measurement is done to one part of the signal amplitude: from 10% to 90%. The steps are listed as following:

- (1) Set the vertical mode to be CH1 or CH2. Feed in the signals to the selected channel.
- (2) Adjust the voltage attenuator and the variable to keep the signal vertical amplitude to be displayed for 5 divisions.
- (3) Adjust the vertical position to make the top and the bottom of the signal be located on the scale of 0% and 100% separately.
- (4) Adjust the sweep speed switch to make the rising edge or the trailing edge be displayed on the screen.
- (5) Adjust the horizontal position to make the 10% of the rising edge be located on a certain vertical scale.
- (6) Measure the horizontal distance between two points from 10% to 90%. If the rising or trailing edge is too fast, the horizontal magnification $\times 10$ can be used to magnify the wave to be 10 times in horizontal direction.
- (7) Use the following formula to work out the rising or trailing time of the wave:
Rising (trailing) time=[Horizontal distance (DIV) \times Sweep Time Factor (Time/Div)]/
Horizontal Magnification Rate

For example: In Fig. 5-8, the distance from 10% to 90% of the rising edge of the wave

is 2.4 divisions, the sweep time factor is 1us/DIV, the horizontal magnification rate is $\times 10$, use the formula to work out the rising time as following:

Rising

$$\text{time}=(1\text{us}/\text{DIV}\times 2.4\text{DIV})/10=0.24\text{us}$$




Fig.5-8

5.3.4 Measurement of Time Difference

Use the following steps to measure the time difference of two relative signals:

- (1) Feed the reference signal and the compared signal into CH1 and CH2 separately.
- (2) Set the vertical mode to be ALT or CHOP according to the frequency.
- (3) Set the trigger source to be the reference signal channel.
- (4) Adjust the voltage attenuator and the variable to display the proper amplitude.
- (5) Adjust the level to get stable wave displayed.
- (6) Adjust TIME/DIV to make sure that there is a horizontal distance between two measuring points of the waves, which is convenient for observation.
- (7) Adjust the vertical position to make the measuring points be located on the horizontal scale in the center of the screen.

Time difference=[Horizontal Distance(DIV) \times Sweep Time Factor (Time/DIV)]/ Horizontal Magnification Rate

For example: In Fig. 5-9, the sweep time factor is 10us/DIV, the horizontal magnification rate is 1, the horizontal distance between two points is 1 division, so: Time Difference=(10us/DIV \times 1DIV)/1=10us.




Fig.5-9




Fig.5-10

5.3.5 Measurement of Phase Difference

You can refer the measurement of time difference for the measurement of phase difference using the following steps:

- (1) Set the controlling knobs as step 1 to 4 mentioned above.
- (2) Adjust the voltage attenuator and the variable to make the displayed amplitude the same.
- (3) Adjust the sweep time switch and the variable to make one cycle of the wave be displayed in 9 divisions. Then one division in the horizontal scale would represent 40 degree ($360^\circ \div 9$).
- (4) Measure the relative horizontal distance between two waves.
- (5) Work out the phase difference between two signals using the following formula: Phase Difference=Horizontal Distance (DIV) \times 40 $^\circ$ /DIV

For example: In Fig. 5-10, the relative horizontal distance between two waves is 1 division, then Phase Difference= 40° /DIV \times 1DIV= 40° .

5.4 TV-V Signal Measuring

The TV-V signals can be displayed on the screen using the oscilloscope. The detailed steps are listed as following:

- (1) Set the vertical mode to be CH1 or CH2. Feed the TV-V signals into the selected channel.
- (2) Set the trigger mode to be TV-V and set the sweep speed switch to be 2ms/DIV.
- (3) Observe whether there is a negative synchronized signal is displayed on the screen. If not, please switch CH1 to CH2 and press CH2 INT in to inverse the positive synchronized TV signal into the negative synchronized TV signal.

(4) Adjust the voltage attenuator and the variable to display proper amplitude.

(5) Horizontal magnification $\times 10$ can be used if necessary.

6. Maintenance

!!! Warning: Only qualified professional can do the maintenance.

6.1 Fuse Replacement

The power indicator would be off and the oscilloscope wouldn't function if the fuse is broken. Usually the fuse would not be open-circuited except there is anything wrong with the circuit. Please check the circuit which would lead to the broken fuse. Replace the broken fuse with the specific one.

!!! Warning: Only the fuse with 250V and corresponding current can be used.

Shut off the oscilloscope from the electric net.

6.2 Cleaning

Please use soft cloth with neutral detergent and water to clean the oscilloscope. Don't spray the detergent directly to the surface of the machine so that the water can be prevented entering inside the machine.

Don't use the detergent which contents chemical matters like gasoline, benzene, toluene(-uol), xylene, acetone and etc.

Don't use the matters like abrasive to clean the machine.

6.3 Standard layout

- | | |
|---------------------|-------------------------------------|
| 1. Operation manual | 1 copy |
| 2. Power cable | 1 pc |
| 3. Probe | 2 pcs |
| 4. Test clamp | 1 pcs (only for specifically model) |

Edu-Lab

Unit 1, Karoo Close
Bexwell Business Park

Bexwell, Norfolk

PE38 9GA

01366 385777 | sales@edulab.com