

MEDISAN 70% HAND SANITIZING GEL

Description

A 70% alcohol hospital specification clear hand sanitizing gel.

Product Features

- Sanitizes the skin without the need to water.
- Conforms to BS EN 14476 enveloped viruses with a 60 second contact time (clean conditions)
- Conforms to BS EN 1276 with a 60 second contact time (dirty conditions)
- Conforms to BS EN 1650 (yeast) with a 60 second contact time (dirty conditions)
- Conforms to BSEN 1500 Chemical Disinfectants and Antiseptics Hygienic Handrub Test Method
- Contains emollients to aid in the prevention of dry skin
- Ideal for use in non-clinical hospital environments, on entering food production areas, or before handling food
- An alcohol-based product for use in any situation where microbial control is important
- Evaporates quickly and is non-tacky

Directions for Use

Apply a small amount to the hands and rub well, paying particular attention to finger nails and the area between the fingers.

Allow hands to air dry.

Important – the use of towels or tissue may lead to recontamination with bacteria.

Application to dry cracked or cut skin should be avoided.

Thoroughly rinse the dispenser with warm water, pumping through the nozzle at least once a month to avoid blockage.

Environmental

- Do not allow large quantities of neat product to enter surface waterways.
- In the event of a large spillage, contact the local water authority or Environment Agency immediately.
- This product uses recyclable packaging. The detergents used in the formulation are biodegradable to EU guidelines under the detergents regulations (EC) No 648/2004.
- It contains no Phosphates.
- None of the components of this product are animal derived.

Storage

Store in cool conditions away from direct sources of heat or sparks. Avoid extremes of temperature Do not allow to freeze. Ensure the cap is tightly fitted after use. Please keep out of reach of children.

Properties

Appearance:

Density @ 20°C:

pH (concentrated):

Shelf Life

30 months if stored as recommended

Available in the following containers

6 x 500m