

# Rust Remover Inhibited 2

## Description

Rust Remover Inhibited 2 is a high-performance rust remover for many industrial applications. It is an ideal product for fast, efficient removal of rust from plant, vehicles, oil tools, scaffolding, etc. It is based on phosphoric acid and contains inhibitors to prevent flash rusting, after use.

Rust Remover Inhibited 2 also protects against biological fouling which is a massive issue for operators globally. Many crustaceans, barnacles and seaweeds can grow on all types of equipment, affecting the efficiency and flow of capacity of operations.

## Application

Immerse corroded items into liquid until all corrosion is removed. Hose off product with fresh water. Depending on severity of the corrosion, it can be used neat, or diluted 4:1 with water. The product is also available as a gel for the safe removal of rust from vertical surfaces and even from overhead locations.

## Dilution

The product can be used as in dip-tank and may be diluted up to 10:1 for this purpose. For severe rusting problems, brush onto surface neat, leave for sufficient contact time for the rust to be removed, then rinse with copious amounts of water.

## Environmental

Dispose of chemicals in accordance with local environmental legislation Local environment registrations.

Rust Remover Inhibited 2 is CEFAS listed with a Gold rating.

## Handling

Normal precautions should be taken for handling chemicals. Protective clothing, eye protection and gloves should be worn. For further, details please refer to the safety data sheet.

## Storage

Store in a cool, dry location between 5°C and 25°C. Protect from freezing. Avoid direct sunlight and extreme temperatures. Store the product sealed in the original container to avoid air exposure. Some products may settle during storage or transportation, ensure thorough mixing before use. See safety datasheet for list of incompatible materials and other relevant information.

## Properties

Appearance: Blue liquid  
Density @ 21°C: 1.19 - 1.21  
pH (concentrated solution): <2

## Shelf Life

Storage: 24 months

## Available in the following containers

